

GLYCANPAC AXH-1 COLUMN

Quick Start

1. Overview

The Thermo Scientific™ GlycanPac AXH-1 column is a high-performance, silica-based HPLC column for simultaneous separation of glycans by charge, size and polarity. It is designed for high-resolution and high-throughput analysis with unique selectivity for biologically important glycans, either labeled or un-labeled, by LC-fluorescence and LC-MS methods.

2. Main features

- Unique glycan selectivity based on charge, size and polarity.
- Excellent resolution for both unlabeled and labeled glycans.
- Useful for both high-resolution glycan profile characterization and easy quantification of glycans based on charge.
- Compatible with fluorescence and MS detection methods.
- High chromatographic efficiency and excellent column stability.

3. Physical data (Table 1)

	GlycanPac AXH-1 Column (3 μm)	GlycanPac AXH-1 Column (1.9 μm)
Column chemistry	WAX and HILIC Mixed-Mode	WAX and HILIC Mixed-Mode
Silica substrate	Spherical, high-purity, porous	Spherical, high-purity, porous
Particle size	3 μm	1.9 μm
Surface area	300 m ² /g	220 m ² /g
Pore size	120 Å	175 Å

4. Specifications and Recommended Operational Parameters (Table 2)

Column Particle size	Column Dimension	P/N	Maximum Pressure (psi)	pH Range	Temperature Limit (°C)	Solvent/Aqueous Compatibility	Recommended Flow Rate (mL/min)	Maximum Flow Rate (mL/min)
1.9 μm	2.1 x 100 mm	082473	7,000	2.0 – 8.0	< 60	Compatible with 0 – 90% aqueous and common HPLC solvents (except acetone)	0.1 - 0.4	0.50
	2.1 x 150 mm	082472	10,000	2.0 – 8.0	< 60		0.1 - 0.4	0.50
	2.1 x 250 mm	082471	15,000	2.0 – 8.0	< 60		0.1 – 0.4	0.50
3 μm	4.6 x 150mm	082468	6,000	2.0 – 8.0	< 60		0.6 – 1.2	1.50
	3,0 x 150 mm	082469	6,000	2.0 – 8.0	< 60		0.3 – 0.6	0.75
	2.1 x 150 mm	082470	6,000	2.0 – 8.0	< 60		0.1 – 0.4	0.50

5. Operational Guidelines

- **All new columns** or any column not in use for longer than 3 days **must** be treated using the cleaning procedure described below in Table 3 before use. Equilibrate the column and perform 3 blank injection prior to analysis of real samples.



NOTE

Buffered solutions must be used for analysis and storage.

- Operate the column within operating specifications (see Table 2 for details).
- Avoid sudden pressure surge.
- Follow the direction of flow is marked on the column.
- Always use guard columns when analyzing samples and replace them before exhausted.
- Column storage: use mobile phase for short-term storage (< 24 hours) and a solution containing 90% acetonitrile and 10% ammonium formate buffer (e.g. 100 mM, pH4.4) for long-term storage (> 24 hours).
- Mobile phase: acetonitrile/ammonium formate buffer (e.g. 100 mM, pH4.4) system is recommended for both LC/fluorescence and LC/MS applications
- Whenever abnormal peak broadening and/or peak tailing is observed, perform the column cleaning described below.




NOTE

Column performance is affected by the contaminants from samples, LC system, and mobile phase. Metal contamination can be witnessed quite often because most LC system and column hardware are made from stainless steel which will bleed out iron over time. When this happens, perform the column cleaning e described below (also in Section 4.7 in the GlycanPac AXH-1 Column Manual) to restore the performance of the column.

6. Cleaning Procedure (Table 3)

Time (min)	Acetonitrile (%)	50 mM sodium pyrophosphate in 100 mM ammonium formate, pH 4.4 (%)	D.I. water (%)	100 mM Ammonium formate buffer, pH 4.4 (%)	Flow (mL/min)		
					2.1 mm dia. column	3.0 mm dia. column	4.6 mm dia. column
0	50	0	0	50	0.25	0.51	1.20
5	20	0	0	80	0.25	0.51	1.20
10	20	80	0	0	0.25	0.51	1.20
35	20	80	0	0	0.25	0.51	1.20
36	20	0	0	80	0.25	0.51	1.20
50	20	0	0	80	0.25	0.51	1.20
60	80	0	0	20	0.25	0.51	1.20
61	78	0	20	2	0.25	0.51	1.20
70	78	0	20	2	0.25	0.51	1.20

These conditions are appropriate for all 100 or 150 mm long columns. Increase the times by 60% for the 250 mm long column.



Performance comparison before and after cleaning of the column for 2-AB labeled glycans from bovine fetuin. Peaks are grouped by sialylation number.

7. Ordering Information (Table 4)

	Particle Size	Column Dimensions	P/N	Required Holder
Analytical	1.9 μm	2.1x250 mm	082521	
		2.1x150 mm	082472	
		2.1x100 mm	082473	
	3 μm	4.6x150mm	082468	
		3.0x150 mm	082469	
Guard pkg. of 2	3 μm	2.1 x 150 mm	082470	
		4.6 x 10 mm	082474	P/N 069580
		3.0 x 10 mm	082475	P/N 069580
		2.1 x 10 mm	082476	P/N 069580