



Prep and Load Platform

ITEX-2 Option

High Sensitivity
Enrichment Technique for
Gas Chromatography



Environmental / Drinking Water
Food / Flavour / Consumer Products
Forensics / Toxicology
Petrochemicals / Polymers
Pharmaceuticals / Residual Solvents

PAL Combi-*xt* General Specifications

System Type

XYZ robot with syringe only concept, no tubing in sample path

Local User Interface

Control panel with 4 function keys, graphical LCD display, unique scroll knob for teach functions

Remote Control

Cycle Composer control software Windows 2000 / XP
Third party instrument drivers for all major GC/GC-MS Systems

Maintenance

Accessibility to all maintenance parts from front
Preventative maintenance kits available

Electrical Control

1x RS232 / 1 x LAN (with optional PAL Upgrade Electronics)
3x TTL Input
2x Opto Coupler Input
2x Relay Output

Power Requirements

100-240V, 120W, 50/60Hz

Environment

4°C - 40°C constant temperature, < 80% humidity (non condensing)

Weight

~ 10kg (without accessories)

Dimension

Length 828mm Depth 385mm Height 575mm

Electrical Safety Standards

CAN/CSA C22.2 No. 61010-1 / ANSI/UL 61010-1 / EN 61010-1

Specifications are subject to change without notice

Sample Capacity*

up to 600	1ml micro vials (78 1ml vials standard)
294	2ml vials (98 2ml vials standard)
96	10ml or 20ml vials
4	deepwell microplates (96/384 wells)
8	standard microplates (96/384 wells)

(* depends on GC model)

GC Mounting Kits

Agilent Technologies 5890 / 6850 / 6890 | 7890
Thermo Scientific GC 8000Top / TRACE GC / Focus GC
Varian GC 3400 / 3600 / 3800 / 3900 / 430 / 450
Shimadzu GC 14 / 17A / 2010 / 2014
Perkin Elmer Autosystem XL / Clarus 400 / Clarus 500 / Clarus 600
GL Sciences GC 353 / 393 / 4000

Order details for ITEX Option (part no. PAL ITEX-2Option)

Description

1pc	ITEX-2 Syringe 1.3mL with M7 x 0.5 Fitting
1pc	Replacement plunger 1.3ml
2pc	ITEX-2 trap TENAX TA 80/100 mesh
1pc	Trap heater incl. electrical connections
1pc	Endplate left side
1pc	Syringe heater side bracket
1pc	CD-ROM including ITEX Cycle (requires Cycle Composer)

Consumables

ITEX-2TrapTXTA	1pc ITEXTrap Tenax TA
ITEX-2TrapTXTA3	Set of 3pcs. ITEXTrap Tenax TA
SYRC ITEX-2.-1.3	1pc replacement ITEX Syringe 1.3mL
PLG ITEX-2.-1.3	Replacement plunger for 1.3mL syringe

Custom filled traps available on request
Please inquire with your local distribution partner

To learn more about the unique PAL-*xt* Series of GC/GC-MS sample injection systems or any of our LC/LC-MS sample handling systems contact your PAL System distributor.

Static Headspace - Liquid Injection - SPME - ITEX Extraction combined in one single instrument

Distributed by:

CTC Analytics

Where design meets performance

CTC Analytics AG
Switzerland
info@ctc.ch
www.palsystem.com

CTC Analytics acknowledges all tradenames and trademarks used as the property of their respective owners





ITEX adsorption step out of a sample vial



ITEX-2 upgrade kit including
2pc. ITEX traps containing proven
industry standard adsorbents

Specifications ITEX-2 Option

Pumping Syringe Size: 1.3ml HD syringe with removable trap

ITEX-2 Trap: Stainless steel material, deactivated by Siltek®:
Needle: Injection Needle gauge 23, Point style 5 (side hole)
Standard Trap Material: 44mg Tenax TA 80/100 mesh

Extraction Speed: selectable from 10µL/s up to 1000µL/s

Extraction Strokes: Selectable from 1 - 999

Extraction Volume: Selectable from 130µL - 1300µL/stroke

Desorption Temperature:
+5°C above ambient - 350°C selectable in 1°C increments

Heating-up rate: up to 12°C/s

Desorption Speed: 1µL/sec. - 500µL/sec.

Pumping Syringe and Trap Cleaning:
Inert gas purging, 30sec. - 3600min.

Heated Pumping Syringe:
+5°C above ambient - 150°C selectable in 1°C increments

Incubator Oven:
6 heated vial positions for 2mL / 10mL / 20mL vials
+5°C above ambient - 200°C selectable in 1°C increments

Agitation:
Interval shaking 250rpm - 750rpm, selectable in 1rpm increments

Incubation Time: Up to 999 minutes selectable in 1 second increments

CTC Analytic's aim is to supply instruments to customers which make the operation of sample processing simple and transparent. In-line with today's lab requirements for productivity, CTC expanded the application range of its GC Injector System PAL Combi-xt introducing the ITEX-2 Option. The ITEX-2 Option consists of an add-on module which can be used with any existing or new PAL Combi-xt System. It performs enrichment of volatile or semi-volatile compounds during headspace analysis. A microtrap filled with adsorbent material, such as Tenax or activated charcoal is placed between the heated PAL Combi-xt Headspace syringe and the syringe needle. Using the HS syringe as a pump, a part of the gaseous phase of the pre-conditioned sample vial is pumped repeatedly through the microtrap. This system setup allows rapid, simple and efficient extraction of volatile and semi-volatile sample compounds. To gain sensitivity simply the number of pumping strokes can be increased or several different vials containing the same sample can be extracted. During thermal desorption into the GC Injector the microtrap is rapidly flash heated and the analytes reach the GC column as a narrow band. No cryofocussing is needed to obtain sharp peaks. To prepare the next extraction, the hot trap is re-conditioned outside the injector with clean purge gas.

ITEX-2 trap material examples

Tenax TA

Volatile and semivolatile compounds, temperature limit of 350°C

Carbotrap/Carbopack

Non-porous graphitized carbon blacks (GCBs)
Hydrophobic properties minimized sample displacement by water

Carbosieve/Carboxen

For very volatile compounds, e.g. Vinylchloride, Freon compounds



ITEX-2 Sample Extraction Procedure

Sample is heated and / or agitated in a sealed sample vial until equilibrium is achieved

The ITEX Trap needle pierces the sealed sample vial and the heated gas tight syringe pumps the gas phase through the ITEX Trap

The loaded ITEX Trap is flash heated up to 350°C and desorbed into the hot GC Injector

After thermal desorption the hot ITEX Trap is cleaned with inert flush gas



ITEX-2 parameter control by Cycle Composer

Method	Inj. Vol.	Trap	First Val.	Last Val.	Count
1. Single GC Injection	2	2w2	1	28	1
2. PAL, Local GC Injection	30	13w13	17	24	3
3. Calibration D	2	2w2	17	24	3
4. PAL, Local GC Injection	30	13w13	17	24	3
5. PAL, Local GC Injection	30	13w13	17	24	3
6. Calibration 6-13	2	2w2	1	16	1
7. Sample#112	2	2w2	1	14	1
8. Sample#113	6	2w2	29	43	1
9. Sample#114	6	2w2	57	78	1
10. Calibration 1-19	4.5	2w2	1	14	4
11. Single GC Injection	4.5	2w2	1	5	4
12. Single GC Injection	4.5	2w2	29	35	4
13. Single GC Injection	4.5	2w2	51	56	4
14. Single GC Injection	4.5	2w2	95	91	4
15. PAL, Local GC Injection	2	2w2	44	55	4

Cycle Composer sample list

Flexible Software Control

Choose between two options to control your PAL Combi-*xt* ITEX-2 Option. For individual application requirements the standalone PC based Windows XP / Vista software Cycle Composer is available.

For single keyboard operation of a whole GC/GC-MS system, the following third party PAL Combi-*xt* drivers are available*.

Vendor


Agilent
Agilent
DataApex
Dionex
Justice Software
Leco
Shimadzu
Thermo Scientific
Varian
Varian
Waters
Waters

Software

ChemStation
EZChrom Elite
Clarity
Chromleon
Chromperfect
ChromaTOF
GCMSsolution
Xcalibur
Star
Galaxie
Masslynx
Empower

* certain drivers may not support the ITEX cycle

EPA 502.2 (Calibration Mix) with ITEX



Comparison of ITEX analysis versus Static Headspace
Sample: Purge and Trap calibration mix
(Restek Cat.No. 30431 502.2 CAL2000 Mega-Mix)

Static Headspace Parameter

60°C / 10min / 1mL sample volume

ITEX Parameter

Extraction Speed: 100µL/sec.

Total Pumping Strokes: 50

Temperature Pumping Syringe / Sample Incubation: 60°C / 10min.

Desorption at 200°C, 15sec. splitless

Chromatography:

Injection: Splitless 15sec. at 250°C / Carrier gas: 0.2bar hydrogen

Column: Rtx-502.2 60m x 0.32mm ID, 1.8µm film

Temperature Program: 40°C - 1min. - 10°C / min to 220°C

Detection: FID 250°C

1 1,1-Dichloroethylene	29 m-Xylene
2 Methylene chloride (dichloromethane)	30 p-Xylene
3 trans-1,2-Dichloroethylene	31 o-Xylene
4 1,1-Dichloroethane	32 Styrene
5 2,2-Dichloropropane	33 Isopropylbenzene
6 cis-1,2-Dichloroethylene	34 Bromoform
7 Chloroform	35 1,1,2,2-Tetrachloroethane
8 Bromochloromethane	36 1,2,3-Trichloropropane
9 1,1,1-Trichloroethane	37 n-Propylbenzene
10 1,1-Dichloropropene	38 Bromobenzene
11 Carbon tetrachloride	39 1,3,5-Trimethylbenzene
12 1,2-Dichloroethane	40 2-Chlorotoluene
13 Benzene	41 4-Chlorotoluene
14 Trichloroethylene	42 tert-Butylbenzene
15 1,2-Dichloropropane	43 1,2,4-Trimethylbenzene
16 Bromodichloromethane	44 sec-Butylbenzene
17 Dibromomethane	45 4-Isopropyltoluene (p-Cymene)
18 cis-1,3-Dichloropropylene	46 1,3-Dichlorobenzene
19 Toluene	47 1,4-Dichlorobenzene
20 trans-1,3-Dichloropropylene	48 n-Butylbenzene
21 1,1,2-Trichloroethane	49 1,2-Dichlorobenzene
22 1,3-Dichloropropane	50 1,2-Dibromo-3-chloropropane
23 Tetrachloroethylene	51 1,2,3-Trichlorobenzene
24 Dibromochloromethane	52 Hexachloro-1,3-butadiene (Hexachlorobutadiene)
25 1,2-Dibromoethane (EDB)	53 Naphthalene
26 Chlorobenzene	54 1,2,3-Trichlorobenzene
27 1,1,1,2-Tetrachloroethane	
28 Ethylbenzene	

Volatiles with ITEX




Figure 1: TIC (m/z 29-400) of Volatile Organic Compounds in Beer. Additional components could be identified due to 10 x higher sensitivity of ITEX compared to Static Headspace.

1 1-Propanol	5 3-methyl-1-butanol	9 3-methyl butyl acetate
2 Ethylacetate	6 2-methyl-1-butanol	10 2-methyl butyl acetate
3 2-methyl-1-propanol	7 2-methyl propyl acetate	
4 Ethyl propanoate	8 Ethyl butyrate	

Beer Ketones with ITEX




Figure 2: Extracted Ion chromatograms for m/z 86 obtained by GC/MS in SIM mode. The Diacetyl Peak can be detected with at least 6 times better S/N value using ITEX rather than Static Headspace. The concentration of diacetyl in this beer sample was in the order of less than 10ppb

Static Headspace Parameter

80°C / 15min / 1ml sample volume

ITEX Parameter

Extraction Speed: 50µL/sec.

Total Pumping Strokes: 10 x 1mL

Temperature Pumping Syringe / Sample Incubation: 80°C / 15min.

Desorption at 250°C

Trap Material: Tenax TA 80/100 mesh

Chromatography:

Injection: Split 1:25 at 250°C

Carrier gas: 200 kPa He at constant pressure

Column: DB-VRX 20m x 0.18mm i.d. / 1µm film

Temperature Program: 40°C - 5min. - 10°C / min to 250°C - 10min.

MSD transfer line: 250°C (17 cm x 110µm i.d. restrictor, 28kPa)

Detection: MS in Scan/SIM Mode

Scan: 29-400 amu

SIM Ions monitored: 43, 57, 86, 100 (50ms dwell time)

- Get P&T sensitivity without the cost of a P&T system
- Rapid & efficient sample enrichment of volatile & semi-volatile compounds in solid, liquid and gaseous samples
- In-tube extraction and direct thermal desorption using proven industry standard adsorbents
- Syringe only concept for transparent sample handling, no sample loops, no transfer lines, no switching valves
- No GC injector modifications, no cryo-focussing necessary
Top mounted on GC's, saves valuable bench space
- Interfaces with any System controlled by all major GC/GC-MS Systems



PAL SYSTEM

www.palsystem.com