

Leica HCS A

High Content Screening Automation
 Technical Documentation

Living up to Life

Leica HCS A

High Content Screening Automation

LAS AF MATRIX M3 Software Platform

Amplify the power of your imaging system by HCS A and convert your confocal into a High Content Screening System. Maximum flexibility is provided for universal applications and easy automation enables efficient screening. Computer Aided Microscopy allows external system control and turns your Leica imaging system into an intelligent microscope.




Leica TCS SP5



Leica TCS SPE



Leica TCS LSI



Leica HCS A

Leica HCS A

Product Overview				
Imaging automation	LAS AF MATRIX Mosaic Advanced			156602501
Remote system control	LAS AF MATRIX Mosaic+Multiwell Advanced			156602502
	LAS AF MATRIX Mosaic Full Version			156602504
	LAS AF MATRIX Multiwell Full Version			156602505
	LAS AF MATRIX Full Version w/o CAM			156602511
	LAS AF MATRIX Developer Entry			156602512
Accessory Tools	LAS AF MATRIX Developer Full w. CAM			156602514
	LAS AF MATRIX Single Object Tracking			156602507
Upgrades	LAS AF MATRIX Pump ²			156602508
	LAS AF MATRIX Z-Drift Compensator			156602509
	LAS AF MATRIX Advanced Trigger ²			156602510
	LAS AF MATRIX Upgrade Mosaic Adv. to Mosaic Full Version			156602515
	LAS AF MATRIX Upgrade Developer Entry to Full Version			156602519
Technical Data				
Platforms	System support ¹		TCS SP5, TCS SPE, TCS LSI	
	LAS AF version		2.2.1 or higher	
Image acquisition	Imaging technologies		True confocal point scanner	
Export formats	Image types		TIF, OME.TIF, LIF	
	Image data		OS platform independent	
Network	Protocol		TCP/IP	
	Administration		Local system admin	
Remote system control	Control via network		Yes, with CAM	
	Control interface		Computer Aided Microscopy, CAM	
Compatible software (examples)	Not included in the packages			
	Programming languages		C++, C#, VB LabVIEW™, MATLAB™	
	Image programming		Adobe Photoshop®	
Limitations	Image analysis software		ImageJ ² with LOCI plugins Definiens® Challenger® MetaMorph®, Leica MM AF®	
	LAS AF MATRIX M3, version 2.2.1		Hardware excluded:	TCS SP5 X, MP, STED, DMI 6000 CFS
			Software excluded:	TCS SMD FCS, FLIM, FLCS FRAP, FRET, Electrophysiology
System Requirements				
Imaging systems Computing	Platform		TCS SP5, TCS SPE, TCS LSI	
	Operating System		Windows® XP, SP2	
	Recommended			
	Processor speed	Mhz	3	
	Memory	GB	4	
	Hard disc	GB	500	
	Network		Yes	
Motorized Stages	Platform	Stand	Part No.	XY-Travel Range³
	TCS SP5	DM	15 6504 150	100 x 80 mm
		DMI	15 6504 141	127 x 83 mm
	TCS SPE	DM	15 6605 200	72 x 38 mm
		DMI	15 6605 201	127 x 83 mm
TCS LSI	LSI-Base	15 6905 202	127 x 83 mm	

Annotations:

¹ Supports technology within the range of the product specification.

² Supported by TCS SP5 only with additional external hardware.

³ Please note that depending on different well plate formats and platforms the outer wells might not be scanned due to mechanical restrictions. Test is recommended.

Open Microscopy Environment (OME) is a multi-site collaborative effort among academic laboratories and a number of commercial entities that produces open tools to support data management for biological light microscopy. Designed to interact with existing commercial software, all OME formats and software are free, and all OME source code is available under GNU public copyleft licenses. OME is developed as a joint project between research-active laboratories at the Dundee, NIA Baltimore and Harvard Medical School and LOCI.

C++ is a programming language standardized by ISO. C# is a programming language developed by Microsoft, Inc. LabVIEW™ is a registered trademark of NI National Instruments Inc.

MATLAB™ is a registered trademark of The MathWorks™, Inc. Adobe Photoshop® is a registered trademark of Adobe Systems® Incorporated.

ImageJ is a public domain Java image processing program inspired by National Institutes of Health, NIH. Windows® XP is a registered trademark of the Microsoft® Corporation.

Definiens® is a registered trademark of Definiens AG. Meta-Morph® is a registered trademark of MDS Analytical Technologies.

“With the user, for the user”

Leica Microsystems

Leica Microsystems operates globally in four divisions, where we rank with the market leaders.

• Life Science Division

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

• Industry Division

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

• Biosystems Division

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

• Medical Division

The Leica Microsystems Medical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

The statement by Ernst Leitz in 1907, “with the user, for the user,” describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Active worldwide

Australia:	North Ryde	Tel. +61 2 8870 3500	Fax +61 2 9878 1055
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Belgium:	Groot Bijgaarden	Tel. +32 2 790 98 50	Fax +32 2 790 98 68
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Ballerup	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Nanterre Cedex	Tel. +33 811 000 664	Fax +33 1 56 05 23 23
Germany:	Wetzlar	Tel. +49 64 41 29 40 00	Fax +49 64 41 29 41 55
Italy:	Milan	Tel. +39 02 574 861	Fax +39 02 574 03392
Japan:	Tokyo	Tel. +81 3 5421 2800	Fax +81 3 5421 2896
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
	Shanghai	Tel. +86 21 6387 6606	Fax +86 21 7387 6698
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Kista	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Heerbrugg	Tel. +41 71 726 34 34	Fax +41 71 726 34 44
United Kingdom:	Milton Keynes	Tel. +44 800 298 2344	Fax +44 1908 246312
USA:	Buffalo Grove/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives in more than 100 countries