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Iron and Magnesium Determination in Meat using Flame Atomic Absorption Spectroscopy

Dr Anastasia Gadzhieva, AA Applications Chemist, Thermo Fisher Scientific, Cambridge, UK

Iron, Copper and Zinc Determination in Wine using Flame Atomic Absorption Spectroscopy


Figure 2: An example of a calibration curve for iron.

The results indicate that these meat samples do not exceed the recommended daily intake (RDI) of 15 mg/kg for iron and 350 mg/kg for magnesium.

In vino veritas, in aqua sanitas.

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Table 1. Instrument settings for the ICE 3300 AAS.

Parameter	Iron	Copper	Zinc
Wavelength	248.3 nm	324.8 nm	213.9 nm
Band pass	0.2 nm	0.5 nm	0.2 nm
Background Correction	D2	D2	D2
Lamp Current	75%	75%	75%
Signal	Continuous	Continuous	Continuous
Flame type	Air-Acetylene	Air-Acetylene	Air-Acetylene
Fuel flow rate	0.9 L/min	1.1 L/min	1.2 L/min
Measurement Time	4s	4s	4s
Replicates	3	3	3


Figure 2: Calibration curve for Iron.

*All data was calculated from 3 replicate readings for each solution, the results are calculated with a view of 1:4 dilution.

Sample ID	Iron (mg/L)	Copper (mg/L)	Zinc (mg/L)
White wine*	1.165	0.027	0.988
Rosé wine*	1.835	0.028	0.514
Red wine*	2.437	0.037	0.575

Cadmium Determination in Crab Meat using Graphite Furnace Atomic Absorption Spectroscopy


Figure 3: The calibration curve for cadmium.

Table 3. Results show the measured concentration of cadmium as well as expected and measured concentrations with percentage spike recovery for three separate crab samples (*all data was calculated from 3 replicate readings for each solution using transient peak height measurements).

Sample	Measured concentration sample (mg/kg)	Expected concentration spike (mg/kg)	Measured concentration spike (mg/kg)	Spike recovery (%)
Tinned crab*	0.0005	0.5005	0.4663	93.2
Dressed crab*	0.9641	1.46	1.4594	100.0
Crab pâté*	0.6143	1.06	1.0919	103.0

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Food Safety Series – Accurate analysis of low levels of mercury in fish by vapor generation AA

The Analysis of Cadmium in Chocolate by Graphite Furnace Atomic Absorption Spectrometry

The Determination of Selenium in Shampoo by Flame Atomic Absorption Spectrometry

Determination of Trace Elements in Rice Products by Flame and Graphite Furnace Atomic Absorption Spectrometry

Rebecca Price, Thermo Fisher Scientific, Cambridge, UK

Application Notes 43019

Method Guide: 40184

Atomic Absorption Full Method Sn in Canned Fruit Juice

The Analysis of Trace Elements in Honey by Flame and Graphite Furnace Atomic Absorption Spectrometry

Rebecca Price, Thermo Fisher Scientific, Cambridge, UK

Application Notes 43060

Key Words

Atomic Absorption, Esse Rice, Toxic Elements.

Key

- Th
- Ab
- cap
- per
- An
- pr
- me
- Th
- all
- fur

Sum

The Spectrometer Flame tool is equipped with a copper atom detector lead. optin accur

Intr

Rice with 600 count between crop essence of health cadm Prefe

Key Words

Atomic Absorption, Flame, Deuterium, Graphite Furnace, Honey, Zeeman

Key Benefits

- The robust flame sample introduction system allows dissolved honey samples to be run without blockage or contamination.
- The advanced furnace autosampler speeds up analysis by automatically preparing the working standards from a single master standard.
- The permanently aligned true dual atomizer enables rapid switching between flame and furnace methods.
- Deuterium and Zeeman background correction offer a flexible solution for the analysis of challenging matrices such as honey.

Summary

The Thermo Scientific iCE 3500 Atomic Absorption Spectrometer is the ideal solution for the analysis of major, minor and toxic elements in honey. The permanently aligned true dual atomizer allows robust and reliable analysis of major elements by flame, followed by accurate and precise determination of minor and toxic elements by graphite furnace.

Introduction

Honey is a sweet and viscous substance produced from the nectar and secretions of plants and flowers. The nectar is transported to a beehive by honey bees, where worker bees then add enzymes to create honey. Most honey is created from a variety of plants and flowers, though in some areas, where a particular plant or flower is in abundance, monofloral honey can be produced, and this is particularly valuable. Honey is typically advertised to the consumer by floral source or geographical location, however many honey products are blended from a variety of sources. This has resulted in a global market with hundreds of types of honey, each with unique taste, color and crystallization properties. In the EU, honey must adhere to strict composition criteria, including sugar, moisture and hydroxymethylfurfural (HMF) content.¹

Analysis of Zinc in Toothpaste using the Thermo Scientific iCE 3000 series AAS

Bhagwesh Sunkar, Dr. G. W. Joshi, Dr. B.A. Patra and Rakesh Jha, Thermo Fisher Scientific, Mumbai, India

Application Note 43159

Key Words

AA, Flame, Teeth, Toothpaste, Zinc

Introduction

Typically our day starts with a mixture of chemicals, toothpaste, a paste or gel dentifrice used with a toothbrush as an accessory to clean and maintain the aesthetics and health of teeth. Toothpaste is used to promote oral hygiene: it serves as an abrasive that aids in removing the dental plaque and food from the teeth, assists in suppressing halitosis (bad breath), and delivers active ingredients such as fluoride or xylitol which help to prevent tooth and gum disease (gingivitis).¹⁾

Although toothpaste was used as long ago as 500 BC in ancient Egypt, China and India, it was not until the 19th century that toothpastes came into general use. Ancient tooth powders used abrasive ingredients such as crushed bone, burned and crushed egg, snail or oyster shells, which were used to clean debris from teeth. First 'modern' toothpastes in the 1800s were homemade, with chalk, soap and salt as common ingredients and were sold in jars either as a powder or paste. A dentist in the US first introduced tooth paste in a tube, in 1892. The most important breakthrough in the history of toothpaste was the addition of fluoride in toothpaste in 1914. The use of fluoride in toothpastes became widespread much later in the 1950s.

Modern toothpaste is a complex formulation that has many more added ingredients since it was first introduced. Extensive work has been carried out to introduce an ingredient which can provide effective protection against various dental conditions.

Triclosan is reported as one of the commonly used important antimicrobial agents used in toothpaste that fight the bacteria which contribute to dental plaque. However, Triclosan induces damage and lesions to the cell wall of bacteria resulting in bacteriolysis (death of the cell) and hence its use is restricted or banned in many countries.²⁾

As an alternative to Triclosan many leading manufacturers prefer to use bacteriostatic agents such as zinc (zinc chloride or zinc citrate) which stop the growth of dental plaque bacteria by inhibiting their metabolism. The combination of bacteriostatic and bactericidal agents as toothpaste ingredients is the most effective one way to fight dental plaque and gum disease.

This application note presents a simple, rugged flame atomic absorption spectrometry method for the analysis of zinc in toothpaste. Random samples of known international & local brands in the market were used for this work.


Figure 3: Correction with QuadLine background correction


Figure 4: Correction with Zeeman background correction

Thermo SCIENTIFIC

The Analysis of Clinical Samples by Atomic Absorption Spectrometry

Krev a sérum

moč

3) Analytes:

The metallic elements measured in clinical analyses fall into three groups according to their role in human metabolism, and they are usually classed as essential, toxic or therapeutic.

Essential elements

- Essential major elements:
Ca, Mg, Na, K
- Essential minor elements:
Zn, Cu, Fe
- Essential trace elements:
Cr, Mn, Mo, Co, V, Se, Ni

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iCAP Q - Modely

Product		iCAP Qa	iCAP Qc	iCAP Qs
				
Sensitivity [kcps/ppb]	⁷ Li	40	50	80
STD Mode	¹¹⁵ In	150	220	400
	²³⁸ U	200	300	500
Sensitivity [kcps/ppb]	⁵⁹ Co	-	30	50
KED Mode				
Oxides [%]	CeO/Ce	< 2	< 2	< 2
Doubly Charged [%]	Ba ⁺⁺ /Ba ⁺	< 3	< 3	< 3
Background m/z 4.5 [cps]	STD	< 1	< 1	< 1
	KED	-	< 0.5	< 0.5

iCAP Q - Modely

Product	iCAP Qa	iCAP Qc	iCAP Qs
			
Short Term Stability (10 min) [%]	< 2	< 2	< 2
Long Term Stability (2 h) [%]	< 3	< 3	< 3
Detection Limits [ppt]			
⁹ Be	< 0.5	< 0.5	< 0.5
¹¹⁵ In	< 0.1	< 0.1	< 0.1
²⁰⁹ Bi	< 0.1	< 0.1	< 0.1
Isotope Ratio Precision [%RSD]			
¹⁰⁷ Ag/ ¹⁰⁹ Ag	< 0.1	< 0.1	< 0.1

Dramatically Different Ion Focusing:

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


RAPID Lens – Simplicity and ion optics done right

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



Mn⁺ Co⁺ As⁺ Se⁺

Mass [u]

Effect of Low Mass Cut-Off on in-cell Interference

ANALYTE	CUT-OFF MASS	POTENTIAL INTERFERENT	PRECURSORS
⁴⁵ Sc	29	¹³ C ¹⁶ O ₂ , ¹² C ¹⁶ O ₂ H, ⁴⁴ CaH, ³² S ¹² CH, ³² S ¹³ C, ³³ S ¹² C	H, C, O, S, Ca
⁴⁷ Ti	32	³¹ P ¹⁶ O, ⁴⁶ CaH, ³⁵ Cl ¹² C, ³² S ¹⁴ NH, ³³ S ¹⁴ N	H, C, N, O, P, S, Cl, Ca
⁴⁹ Ti	33	³¹ P ¹⁸ O, ⁴⁸ CaH, ³⁵ Cl ¹⁴ N, ³⁷ Cl ¹² C, ³² S ¹⁶ OH, ³³ S ¹⁶ O	H, C, N, O, P, S, Cl, Ca
⁵⁰ Ti	34	³⁴ S ¹⁶ O, ³² S ¹⁸ O, ³⁵ Cl ¹⁴ NH, ³⁷ Cl ¹² CH	H, C, N, O, S, Cl
⁵¹ V	35	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, O, N, S, Cl
⁵² Cr	36	³⁶ Ar ¹⁶ O, ⁴⁰ Ar ¹² C, ³⁵ Cl ¹⁶ OH, ³⁷ Cl ¹⁴ NH, ³⁴ S ¹⁸ O	H, C, O, N, S, Cl, Ar
⁵⁵ Mn	39	³⁷ Cl ¹⁸ O, ²³ Na ³² S, ²³ Na ³¹ PH	H, O, Na, P, S, Cl, Ar
⁵⁶ Fe	39	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca
⁵⁷ Fe	40	⁴⁰ Ar ¹⁶ OH, ⁴⁰ Ca ¹⁶ OH	H, O, Ar, Ca
⁵⁸ Ni	41	⁴⁰ Ar ¹⁸ O, ⁴⁰ Ca ¹⁸ O, ²³ Na ³⁵ Cl	O, Na, Cl, Ar, Ca
⁵⁹ Co	42	⁴⁰ Ar ¹⁸ OH, ⁴³ Ca ¹⁶ O, ²³ Na ³⁵ ClH	H, O, Na, Cl, Ar, Ca
⁶⁰ Ni	43	⁴⁴ Ca ¹⁶ O, ²³ Na ³⁷ Cl	O, Na, Cl, Ca
⁶¹ Ni	44	⁴⁴ Ca ¹⁶ OH, ³⁸ Ar ²³ Na, ²³ Na ³⁷ ClH	H, O, Na, Cl, Ca
⁶³ Cu	45	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ¹² C ¹⁴ N ³⁷ Cl, ³¹ P ³² S, ³¹ P ¹⁶ O ₂	C, N, O, Na, P, S, Cl
⁶⁴ Zn	46	³² S ¹⁶ O ₂ , ³² S ₂ , ³⁶ Ar ¹² C ¹⁶ O, ³⁸ Ar ¹² C ¹⁴ N, ⁴⁸ Ca ¹⁶ O	C, N, O, S, Ar, Ca
⁶⁵ Cu	47	³² S ¹⁶ O ₂ H, ³² S ₂ H, ¹⁴ N ¹⁶ O ³⁵ Cl, ⁴⁸ Ca ¹⁶ OH	H, N, O, S, Cl, Ca
⁶⁶ Zn	47	³⁴ S ¹⁶ O, ³² S ³⁴ S, ³³ S, ⁴⁸ C, ¹⁸ O	O, C, S
⁶⁹ Ga	47	³² S ¹⁸ O ₂ H, ³⁴ S ₂ H, ³⁷ Cl ¹⁶ O ₂	H, O, S, Cl
⁷⁰ Zn	47	³⁴ S ¹⁸ O ₂ , ³⁵ Cl ₂	O, S, Cl
⁷⁵ As	47	⁴⁰ Ar ³⁴ SH, ⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ³⁷ Cl ₂ H	H, S, Cl, Ca, Ar
⁷⁷ Se	47	⁴⁰ Ar ³⁷ Cl, ⁴⁰ Ca ³⁷ Cl	Cl, Ca, Ar
⁷⁸ Se	47	⁴⁰ Ar ³⁸ Ar	Ar
⁸⁰ Se	47	⁴⁰ Ar ₂ , ⁴⁰ Ca ₂ , ⁴⁰ Ar ⁴⁰ Ca, ³² S ₂ ¹⁶ O, ³² S ¹⁶ O ₃	O, S, Ar, Ca


Effect of Low Mass Cut-Off on in-cell Interference Formation

ANALYTE	CUT-OFF-MASS	POTENTIAL INTERFERENT	PRECURSORS
⁴⁵ Sc	29	¹³ C ¹⁶ O ₂ , ¹² C ¹⁶ O ₂ H, ⁴⁴ CaH, ³² S ¹² CH, ³² S ¹³ C, ³³ S ¹² C	H, C, O, S, Ca
⁴⁷ Ti	32	³¹ P ¹⁶ O, ⁴⁶ CaH, ³⁵ Cl ¹² C, ³² S ¹⁴ NH, ³³ S ¹⁴ N	H, C, N, O, P, S, Cl, Ca
⁴⁹ Ti	33	³¹ P ¹⁸ O, ⁴⁸ CaH, ³⁵ Cl ¹⁴ N, ³⁷ Cl ¹² C, ³² S ¹⁶ OH, ³³ S ¹⁶ O	H, C, N, O, P, S, Cl, Ca
⁵⁰ Ti	34	³⁴ S ¹⁶ O, ³² S ¹⁸ O, ³⁵ Cl ¹⁴ NH, ³⁷ Cl ¹² CH	H, C, N, O, S, Cl
⁵¹ V	35	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, O, N, S, Cl
⁵² Cr	36	³⁶ Ar ¹⁶ O, ⁴⁰ Ar ¹² C, ³⁵ Cl ¹⁶ OH, ³⁷ Cl ¹⁴ NH, ³⁴ S ¹⁸ O	H, C, O, N, S, Cl, Ar
⁵⁵ Mn	39	³⁷ Cl ¹⁸ O, ²³ Na ³² S, ²³ Na ³¹ PH	H, O, Na, P, S, Cl, Ar
⁵⁶ Fe	39	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca
⁵⁷ Fe	40	⁴⁰ Ar ¹⁶ OH, ⁴⁰ Ca ¹⁶ OH	H, O, Ar, Ca
⁵⁸ Ni	41	⁴⁰ Ar ¹⁸ O, ⁴⁰ Ca ¹⁸ O, ²³ Na ³⁵ Cl	O, Na, Cl, Ar, Ca
⁵⁹ Co	42	⁴⁰ Ar ¹⁸ OH, ⁴³ Ca ¹⁶ O, ²³ Na ³⁵ ClH	H, O, Na, Cl, Ar, Ca
⁶⁰ Ni	43	⁴⁴ Ca ¹⁶ O, ²³ Na ³⁷ Cl	O, Na, Cl, Ca
⁶¹ Ni	44	⁴⁴ Ca ¹⁶ OH, ³⁸ Ar ²³ Na, ²³ Na ³⁷ ClH	H, O, Na, Cl, Ca
⁶³ Cu	45	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ¹² C ¹⁴ N ³⁷ Cl, ³¹ P ³² S, ³¹ P ¹⁶ O ₂	C, N, O, Na, P, S, Cl
⁶⁴ Zn	46	³² S ¹⁶ O ₂ , ³² S ₂ , ³⁶ Ar ¹² C ¹⁶ O, ³⁸ Ar ¹² C ¹⁴ N, ⁴⁸ Ca ¹⁶ O	C, N, O, S, Ar, Ca
⁶⁵ Cu	47	³² S ¹⁶ O ₂ H, ³² S ₂ H, ¹⁴ N ¹⁶ O ³⁵ Cl, ⁴⁸ Ca ¹⁶ OH	H, N, O, S, Cl, Ca
⁶⁶ Zn	47	³⁴ S ¹⁶ O, ³² S ³⁴ S, ³³ S, ⁴⁸ C, ¹⁸ O	O, C, S
⁶⁷ Zn	47	³² S ³⁴ SH, ³³ S ₂ H, ⁴⁸ Ca ¹⁸ OH, ¹⁴ N ¹⁶ O ³⁷ Cl, ³⁵ Cl ¹⁶ O ₂	H, N, O, S, Cl, Ca
⁶⁸ Zn	47	³² S ¹⁸ O ₂ , ³⁴ S ₂	O, S
⁶⁹ Ga	47	³² S ¹⁸ O ₂ H, ³⁴ S ₂ H, ³⁷ Cl ¹⁶ O ₂	H, O, S, Cl
⁷⁰ Zn	47	³⁴ S ¹⁸ O ₂ , ³⁵ Cl ₂	O, S, Cl
⁷⁵ As	47	⁴⁰ Ar ³⁴ SH, ⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ³⁷ Cl ₂ H	H, S, Cl, Ca, Ar
⁷⁷ Se	47	⁴⁰ Ar ³⁷ Cl, ⁴⁰ Ca ³⁷ Cl	Cl, Ca, Ar
⁷⁸ Se	47	⁴⁰ Ar ³⁸ Ar	Ar
⁸⁰ Se	47	⁴⁰ Ar ₂ , ⁴⁰ Ca ₂ , ⁴⁰ Ar ⁴⁰ Ca, ³² S ₂ ¹⁶ O, ³² S ¹⁶ O ₃	O, S, Ar, Ca

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BEC = Background Equivalent Concentration

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Vysoká citlivost

Vysoká tolerance vůči matici


- UNIQUE Cone Inserts
 - Allows variation of central channel length
 - Interchangeable on iCAP Qc to vary instrument properties

Insert	Channel length	Properties
High Matrix	3.5 mm	High matrix tolerance
High Sensitivity	2.8 mm	High sensitivity

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Lets look at the sensitivity specification and mass response for the iCAP Qc with High Matrix insert:


**Sensitivity
Performance
High Sensitivity**


$$= 1.5 \times$$

**Sensitivity
Performance
High Matrix**


- Now add the sensitivity performance and mass response for the iCAP Qc with High Sensitivity insert:

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

The iCAP Qc typical performance with the High Sensitivity interface

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







Trace(s) of Analysis No 4 : "Level 4"


- Plasma stable in reverse – phase LC applications with ACN gradients as mobile phase with O₂ additional plasma gas
- Software:
 - Qtegra plug-in in place for LC control and triggering

LC spojení (Dionex)

Alkylcíny

Trace(s) of Analysis No 6 : "5 ppb"

- Dionex ICS-5000+ analytical quaternary pump
 - Reverse phase LC: 10% acetic acid with 10-90% ACN gradient mobile phase @ 200 $\mu\text{l}\cdot\text{min}^{-1}$

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Aplikační listy iCAP Q + IC Dionex

Total and speciation analysis of Mercury in contact lens solutions by ICP-MS

Daniel Kütischer, Julian Wills, Lother Rottmann, Thermo Fisher Scientific, Germany

Key Words
ICAP Q, ICP-MS, ICS-5000, ion chromatography, ThermoAnal, Hg, speciation

Goal
To develop a fully quantitative method for the determination of both total Hg and Hg species in contact lens solutions.

Introduction
While there is continual awareness regarding exposure to mercury (Hg) sources in general and MeHg in particular due to its presence in food samples such as fish, less interest is paid to the potential risk from ethylmercury (EtHg, or EtHgX). One of the main reasons for this is the faster degradation and consequently excretion of EtHg in the human body that results in considerably lower chronic toxicity. There remain however potential sources where acute intake of EtHg can occur, for example as a consequence of exposure to disinfectant. Thimerosal is used as a bactericide in multi-dose (typical concentration 0.001 to 0.01%) and in other health related products such as eye drops or contact lens solutions. The compound hydrolyses in aqueous solution to form EtHg and thioacetate which is an effective bactericide. Although no direct correlation between thimerosal usage and potential health risks has been established, the use of thimerosal has been reduced in both Germany and the USA. In use it is still permitted in multi dose vaccines and contact lens solutions at concentrations of up to 100 and 70 mg/kg respectively.

For total Hg determination, gold was added to the samples and the rinse solution to minimize memory effects from the sample introduction system.

Sample and Calibration Solution Preparation
Three different commercially available contact lens solutions were prepared for total Hg and Hg speciation analysis. For total mercury analysis, the contact lens solutions were analyzed after a 2000-fold dilution in 2% HNO₃ / 0.5% HCl. However, as no detectable signal for mercury was found, a lower dilution factor (20-fold) was employed. For Hg speciation analysis, calibration standards were prepared in a matrix solution containing 0.5% NaCl and 0.01% EDTA to mimic the matrix of the contact lens solutions and to promote the formation of the same mercury complexes as in the sample matrix (such as [HgCl₂]²⁺). The standards and the contact lens solutions were then diluted 2000-fold with ultra high purity water prior to injection.

Thermo SCIENTIFIC

IC-ICP-MS speciation analysis of As in Organic Brown Rice Syrup (OBRS) using the Thermo Scientific iCAP Q ICP-MS

Daniel Kütischer, Julian Wills and Lother Rottmann, Thermo Fisher Scientific, Germany and Delfel Jansen, Thermo Fisher Scientific, Switzerland

Key Words
ICAP Q, ICP-MS, ICS-5000, ion chromatography, As, speciation, food

Goal
The analysis of organic brown rice samples (OBRS) for their total arsenic content by ICP-MS and their species concentrations by IC-ICP-MS.

Introduction
Media reports and scientific publications on the determination of arsenic (As) in foodstuffs have sparked renewed interest from consumer groups and politicians leading to responses from national regulatory bodies. Following the publication of a report on high As levels in organic brown rice syrup (OBRS, an ingredient in a variety of organic foods) the United States Food and Drug Administration (FDA) stated that it was carrying out a study on As in rice and rice products that is due to report later in 2012.

Based on previous work¹, OBRS samples were analyzed for their total arsenic content by ICP-MS and then subsequently by IC-ICP-MS to determine the concentration of six arsenic species: the two toxic inorganic species (As(III) and As(V)), and four organic species that are considered to be harmless.

Instrument configuration
All measurements – both total As and speciation – were carried out using the Thermo Scientific iCAP Q ICP-MS. All As analyses were performed in He-KED (Kinetic Energy Discrimination) mode, efficiently reduce polyatomic interferences of (monoisotopic) arsenic at m/z 77. Separations were carried out on Dionex ICS-5000 using the iCAP Q sensitive.

Sample Preparation
Three different OBRS samples were sourced and analyzed. For the total As analysis, the dried wheat digestion method described by Jackson et al. was used. Preparation of the OBRS samples for As analysis was achieved by taking 1.5 g of OBRS, ml of 0.28 M HNO₃ and refluxing for 30 min.

Total Arsenic Results
A calibration curve for the He-KED analysis of shown in Figure 1.

Coupling of an Inert Ion Chromatographic System with ICP-Q-MS for Robust and Accurate Elemental Speciation

Daniel Kütischer, Lother Rottmann, Julian Wills, Shona McSherry, Thermo Fisher Scientific, Germany and Delfel Jansen, Thermo Fisher Scientific, Switzerland

Executive Summary
ICP-Q-MS analytical systems for determination of arsenic and selenium species (As and Se) in complex matrices such as biological samples require the use of ion chromatography (IC) for speciation. The Thermo Scientific iCAP Q ICP-MS is a highly sensitive and accurate elemental speciation system. The Thermo Scientific Dionex ICS-5000 ion chromatography system is a highly sensitive and accurate elemental speciation system.

Overview
Previous to this work, the coupling of the Thermo ICS-5000 ion chromatography system to the Thermo iCAP Q ICP-MS was not possible. This work demonstrates the successful coupling of the Thermo ICS-5000 ion chromatography system to the Thermo iCAP Q ICP-MS. The Thermo ICS-5000 ion chromatography system is a highly sensitive and accurate elemental speciation system. The Thermo iCAP Q ICP-MS is a highly sensitive and accurate elemental speciation system.

Key Words
ICAP Q, ICP-MS, ICS-5000, ICP-MS, ion chromatography, As, Se, speciation

Introduction
The Thermo Scientific iCAP Q ICP-MS is a highly sensitive and accurate elemental speciation system. The Thermo Scientific Dionex ICS-5000 ion chromatography system is a highly sensitive and accurate elemental speciation system.

IC-ICP-MS speciation analysis of As in apple juice using the Thermo Scientific iCAP Q ICP-MS

Daniel Kütischer, Shona McSherry, Julian Wills, Thermo Fisher Scientific, Germany, Delfel Jansen, Thermo Fisher Scientific, Switzerland

Key Words
ICAP Q, Apple Juice, As, speciation, ion chromatography, ICS-5000

Goal
Apple juice samples are analyzed by IC-ICP-MS to determine the concentration of six arsenic species: two inorganic and toxic species (As (III) and As (V)) and four organic species (arsenobetaine (AsB), arsenocholine (AsC), monomethylarsinic acid (MMA) and dimethylarsinic acid (DMA)).

Introduction
Interest in the determination of arsenic (As) species in fruit juices has been triggered by media reports in the US claiming that some apple juices may contain high amounts of arsenic¹. Since fruit juices are frequently consumed, especially by young children, the reported values led to serious concerns about the safety of such products. However, in the experiments conducted, only the total arsenic concentration was assessed, but no detailed investigation of the chemical form of the element was carried out. This is an important distinction since the inorganic forms of arsenic (As (III) and As (V)) are highly toxic, while the organic forms (e.g. arsenobetaine) are not considered to be toxic. Typical levels of total arsenic found in apple juice are lower than the US EPA drinking maximum contaminant level (MCL) of 10 µg/L in apple juice is generally considered safe and is currently not regulated. However, as a consequence, the FDA is currently reviewing data, and may eventually tighten its current guidelines.

Instrument o Chromatograph Thermo Scientific system. Due to a this system is pe- studies. A Thermo column (2 mm I this study, as it a arsenic species' separation of the Thermo Scientific sensitive and ad- proprietary 99' provides both h background wh use of a He pres Discrimination) interferences in i KED mode also lower background an l analysis of (non possibility to use therefore the fl reduction QCd minimize polyat would otherwise

General ana The iCAP Q IC cyclic HFA sp (Elemental Scen nebulizer has a with I.C fittings analytes. A dem I.D. quartz inject for optimum per supplied water


Speciation analysis of Cr (III) and Cr (VI) in drinking waters using anion exchange chromatography coupled to the Thermo Scientific iCAP Q ICP-MS

Daniel Kütischer, Shona McSherry, Julian Wills, Thermo Fisher Scientific, Germany, Delfel Jansen, Thermo Fisher Scientific, Switzerland

Key Words
ICAP Q, Cr, speciation, ion chromatography, Drinking water, ICS-5000

Goal
To develop a sensitive, robust and high throughput method for the trace level analysis of Cr (III) and Cr (VI) species in natural waters using IC-ICP-MS.

Introduction
Due to its widespread use in industrial applications such as chromium plating, dye manufacturing and preservation of wood and leather materials, chromium concentrations in environmental samples are monitored on a routine basis. Both the United States EPA and the European Union have specified maximum admissible chromium concentrations in their respective drinking water directives. As with many other trace elements, chromium (Cr) is typically found in more than one chemical form, each of which with different chemical properties and behavior, such as bioavailability and toxicity. For chromium, Cr (III) is essential to human beings and involved in different processes in the body while Cr (VI) is highly toxic. Total Cr content therefore, in, for example, a drinking water sample does not provide sufficient information to evaluate potential hazards to populations exposed to it. In order to provide this critical information a supporting speciation analysis is required to determine the amounts of the different Cr species in the sample. The speciation analysis of Cr however is a challenging task, since the stability of different Cr species is easily affected by conditions during sample collection and treatment¹. For example, low pH values may lead to the degradation of Cr (VI) to Cr (III) due to the increased redox potential, while high pH values may lead to the precipitation of Cr (III) as Cr(OH)₃. An additional difficulty in the accurate speciation analysis of Cr by ICP-MS are the numerous spectral interferences (e.g. ²³⁹Po/⁹¹P to ⁵²Cr/⁵²Cr) on the most abundant chromium isotope, ⁵²Cr.

Thermo SCIENTIFIC

Thermo SCIENTIFIC

Application Note 1000

60557-0000-0000-0000


IC-ICP-MS for the Speciation of Cr (III) and Cr (VI)

Column	Dionex AG-7 (2x50mm)
---------------	-----------------------------

Mobile Phase	0.35 mol L ⁻¹ HNO ₃
--------------	---



Flow rate	0.4 mL min ⁻¹
-----------	--------------------------

Injection volume	20 μL
------------------	-------


IC-ICP-MS for the Speciation of Cr (III) and Cr (VI)

- Both chromium isotopes, $^{52}, ^{53}\text{Cr}$, can be monitored
- Cr(III) and Cr(VI) are completely baseline separated
- Limits of detection:
 - **0.20 ng L⁻¹ [Cr (VI)]**
 - **0.38 ng L⁻¹ [Cr (III)]**
- Isocratic chromatography performed in less than 3 minutes

Speciation of Cr in Drinking Water

***Analysis of locally sourced tap water (Bremen):
Only Cr (VI) was detected, Concentration 42.5 ng L⁻¹
Injection without further sample treatment***


Chromium Speciation Analysis of Drinking Water

- A sensitive and fast method based on ion chromatography for the speciation of Cr has been developed
 - No need for the addition of EDTA or other charge pairing agents necessary, no additional sample preparation, no additional spectral interferences in the mass spectra
 - The IC eluent, dilute HNO_3 , is relatively cheap, easily sourced and allows detection of Cr species in the single digit ppt range
 - Reduced column dimensions allow reduction of the mobile phase flow rate to 0.4 mL min^{-1}

LC Coupling (Dionex) – As Species in Apple Juice




IC-ICP-MS for As Speciation

Column	Dionex AS-7 (2x250mm)
Mobile Phase	Gradient elution: 20 to 200 mM ammonium carbonate
Flow rate	0.3 mL min ⁻¹
Injection volume	20 µL


IC-ICP-MS for As Speciation

- 6 species detected
- ~7000 cps / ppb
- ~15 minute analysis
- Narrow peak shape improves S/N ratio (10 to 20s)
- Reproducibility < 4%
- Reduced solvent consumption without comprising separation efficiency

Arsenic Speciation Analysis of Apple Juice

- Low concentration of the individual species after dilution

	DMA [ng g ⁻¹]	As (III) [ng g ⁻¹]	MMA [ng g ⁻¹]	As(V) [ng g ⁻¹]	Sum [ng g ⁻¹]	Total As [ng g ⁻¹]
Juice 3	-	0.5 ± 0.01	-	0.8 ± 0.01	1.3	1.7 ± 0.05
Juice 4	0.4 ± 0.05	0.3 ± 0.01	0.1 ± 0.05	0.7 ± 0.01	1.5	1.8 ± 0.05

- Method detection limits (MDL): 0.001 ng g⁻¹ As(V) to 0.011 ng g⁻¹ (MMA)

ICP-OES (iCAP 7000)

NEW

low cost
ICP-OES analysis

The Concept going forward

iCAP 6500

**Fastest/Best
Productivity, and
flexibility through S/W
plug-ins**

iCAP 6300

**Faster/More Productive
for Routine QA/QC and
easier to use**

iCAP 6200

**Entry level ICP –
easier to use**

**Qtegra -
New
Common
Software
Platform
-
Easier to
use and
more plug-
ins for
advanced
flexibility**

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
- iCAP 7400
Radial nebo Dual


- iCAP 7600
Radial nebo Dual

Sprint valve system – How does it work?


Sprint valve system – How does it work?


Inject

The iCAP 7000 Series ICP-OES Core Technologies

Optical design

EMT torch & Duo viewing

Thermo
SCIENTIFIC

iCAP 7000 SERIES
ICP Spectrometer


Clip-in sample intro systems

Drain sensor

CID detector

EMT Torch for Ease of Use for Routine Maintenance

TORCH ORIENTATION LOCK

- Auto alignment of the torch in the torch box
- Automatically establishes robust plasma gas connections

SCREW-THREADED CENTRE TUBE HOLDER

- Simple routine maintenance operations with the plasma on!

Still the best torch design and common ICP-OES/ICP-MS


Ceramic *D-Torch* Accessories – Key Features

- Fully demountable design for cost effective replacement of parts
- Identical geometry to EMT (parallel tube design for optimized gas flows and sample channel stability)
- Interchangeable with EMT design in any iCAP 7000 Series ICP
- Silicon nitride (Sialon) ceramic outer tube – highly durable material
- Alumina intermediate (auxiliary) tube – excellent chemical and temperature resistance
- Compatible with EMT centre tubes to enable maximum application flexibility



Dlouhodobá stabilita iCAP 7400




- Offers flexibility with the ability to have greater control of gas flows for method development of complex samples
 - Ideal for high salt samples
 - Ideal for organic samples
- Offers the benefits of MFC in terms of better long term stability
 - Ideal for long analytical runs

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