




pragolab



VIALKY A PRÍSLUŠENSTVO

Drieňová 34/1712
821 02 Bratislava

Tel.: +421 243 294 436
bratislava@pragolab.sk

www.pragolab.sk



- Skrutkovacie vialky
- Inzerty
- Krimpovacie vialky
- Headspace
- EPA vialky
- Skladovanie
- Certifikované kity
- Striekačkové filtre
- Striekačky

BRATISLAVA | Pragolab s.r.o. Tel.: +421 243 294 436
Drieňová 34/1712 Email: bratislava@pragolab.sk
821 02 Bratislava

PRAHA | Pragolab s.r.o. Tel.: +420 284 813 020
Nad Krocinkou 55/285 Email: pragolab@pragolab.cz
190 00 Praha 9 - Prosek

BRNO | Pragolab s.r.o. Tel.: +420 515 908 772
Jamborova 32/3181 Email: bmo@pragolab.cz
615 00 Brno - Židenice

Short Thread Vials ND9



PGL190273

1.5 ml Short Thread Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, label ad filling lines

Pack of 1000 pcs.



PGL190274

1.5ml Short Thread Vial, 32 x 11.6mm, amber glass, 1st hydrolytic class, wide opening, label and filling lines

Pack of 1000 pcs.



PGL190302

1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, label and filling lines, with overwind-barrier


Pack of 1000 pcs.



PGL190182

9mm Combination Seal: PP Short Thread Cap, blue, with centre hole; RedRubber / PTFE beige, 45° shore A, 1.0mm

Pack of 1000 pcs.



PGL190155

UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1.0mm

Pack of 1000 pcs.



PGL190159

9mm Combination Seal: PP Short Thread Cap, blue, centre hole; Silicone white/PTFE blue, 55° shore A, 1.0mm, slit

Pack of 1000 pcs.



PGL190682

9mm Combination Seal, bonded: PP Short Thread Cap, blue, centre hole; Silicone white/PTFE beige, 45° shore A, 1.0mm

Pack of 1000 pcs.



PGL190683

9mm Combination Seal, bonded: PP Short Thread Cap, blue, centre hole; Silicone white/PTFE beige, 45° shore A, 1.0mm, slit

Pack of 1000 pcs.

Microsampling Vials ND9



PGL190298

Short Thread Vial with integrated Micro-Insert, 32 x 11.6mm, clear glass, 1st hydrolytic class, Base bonded

Pack of 1000 pcs.



PGL190299

Short Thread Vial with integrated Micro-Insert, 32 x 11.6mm, amber glass, 1st hydrolytic class, Base bonded

Pack of 1000 pcs.



PGL190276

1.1ml Microliter Short Thread Vial ND9, 32 x 11.6mm, clear glass, 1st hydrolytic class

Pack of 1000 pcs.



PGL190370

1.5ml PP Short Thread Vial, 32 x 11.6mm, transparent, with filling lines

Pack of 1000 pcs.



PGL190374

0.7ml PP Short Thread Micro-Vial, 32 x 11.6mm, transparent

Pack of 1000 pcs.



PGL190366

0.3ml PP Short Thread Micro-Vial, 32 x 11.6mm, transparent

Pack of 1000 pcs.

Micro-Inserts for Short Thread Vials ND9 with wide opening

**PGL190020**

0.1ml Micro-Insert, 31 x 6mm, clear glass, 1st hydrolytic class, 15mm top

Pack of 1000 pcs.

**PGL190022**

0.1ml Micro-Insert, 31 x 6mm, clear glass, 1st hydrolytic class, 12mm top

Pack of 1000 pcs.

**PGL190023**

0.1ml Micro-Insert, 29 x 5.7mm, clear glass, with attached Plastic Spring

Pack of 1000 pcs.

**PGL190024**

0.2ml Micro-Insert, 31 x 6mm, clear glass, 1st hydrolytic class, flat bottom

Pack of 1000 pcs.

**PGL190026**

0.1ml PP Micro-Insert, 29 x 6mm, clear, 10mm top, with filling lines

Pack of 1000 pcs.

**PGL190027**

0.1ml PP Micro-Insert, 29 x 6mm, clear, 10mm top, filling lines and attached Plastic Spring

Pack of 1000 pcs.

**PGL190028**

0.2ml PP Micro-Insert, 31 x 6mm, clear, flat bottom

Pack of 1000 pcs.

Crimp Neck Vials ND11



PGL190269

1.5ml Crimp Neck Vial,
32 x 11.6mm, clear glass, 1st
hydrolytic class, wide open-
ing, label and filling lines

Pack of 1000 pcs.



PGL190270

1.5ml Crimp Neck Vial,
32 x 11.6mm, amber glass,
1st hydrolytic class, wide
opening, label and filling lines

Pack of 1000 pcs.



PGL190227

11mm Combination Seal:
Aluminium Cap, clear
lacquered, centre hole;
Natural Rubber red-orange/
TEF transparent, 60° shore A,
1.0mm / 100 per bag

Pack of 1000 pcs.



PGL190246

11mm Combination Seal:
Aluminium Cap, clear
lacquered, with centre hole;
Red Rubber / PTFE beige, 45°
shore A, 1.0mm

Pack of 1000 pcs.



PGL190228

UltraClean Closure: 11mm
Aluminium Cap, clear lac-
quered, centre hole; Silicone
white/PTFE red, 45° shore
A, 1.3mm

Pack of 1000 pcs.

Crimping and Decapping Tools



PGL190253

11mm Crimper

Pack of 1 pc.



PGL190257

11mm Decapper

Pack of 1 pc.



PGL190213

11mm Electronic Crimper,
110 - 240V

Pack of 1 pc.



PGL190214

11mm Electronic Decapper,
110 - 240V

Pack of 1 pc.


Crimp Neck Headspace-Vials and Aluminum Crimp Seals ND20



PGL190593

20ml Headspace-Vial, 75.5 x 22.5mm, clear glass, 1st hydrolytic class, DIN Crimp Neck, long neck, rounded bottom

Pack of 1000 pcs.



PGL190542

UltraClean Closure: 20mm Aluminium Cap, plain, centre hole; Silicone transparent blue/PTFE white, 45° shore A, 3.0mm

Pack of 1000 pcs.



PGL190559

20mm Combination Seal: Aluminium Crimp Cap, plain, centre hole; Silicone white/PTFE beige, 45° shore A, 3.2mm

Pack of 1000 pcs



PGL190566

UltraClean Closure: 20mm Magnetic Bimetal Cap, red lacquered, 8mm centre hole; Silicone transparent blue/PTFE transparent, 45° shore A, 3.0mm

Pack of 1000 pcs.



PGL190513

20ml Precision Thread Headspace-Vial, 75.5 x 22.5mm, clear glass, 1st hydrolytic class, rounded bottom (for MAGNETIC screw caps)

Pack of 1000 pcs.



PGL190496

UltraClean Closure: 18mm Magnetic Universal Screw Cap, silver, centre hole; Silicone transparent blue/PTFE white, 45° shore A, 1.3mm

Pack of 1000 pcs.

Crimping and Decapping Tools



PGL190572

20mm Crimper

Pack of 1 pc.



PGL190579

20mm Decapper

Pack of 1 pc.



PGL190525

20mm Electronic Crimper, 110 - 240V

Pack of 1 pc.



PGL190527

20mm Electronic Decapper, 110 - 240V

Pack of 1 pc.

(EPA) Screw Neck Vials ND24



PGL190628

20ml EPA Screw Neck Vial,
57 x 27.5mm, clear glass, 1st
hydrolytic class



PGL190627

40ml EPA Screw Neck Vial,
95 x 27.5mm, clear glass, 1st
hydrolytic class



PGL190643

24mm Combination Seal:
PP Screw Cap, white, centre
hole; Silicone white/PTFE
beige, 45° shore A, 3.2mm,
EPA-Quality




PGL190645

24mm Combination Seal:
PP Screw Cap, white, closed
top; Butyl red/PTFE grey, 55°
shore A, 2.5mm



PGL190612

22mm Septa, Butyl red/PTFE
grey, 55° shore A, 2.5mm



PGL190609

22mm Septa, Silicone white/
PTFE blue, 55° shore A,
1.5mm, cross-slitted

Pack of 1000 pcs.

Pack of 1000 pcs.


Pack of 1000 pcs.

Pack of 1000 pcs.

Pack of 1000 pcs.


Pack of 1000 pcs.

PP Storage Boxes for 1,5 ml Sample Vials




PGL190401

Blue



PGL190402

Orange



PGL190403

Pink



PGL190404

Yellow



PGL190405

Green



PGL190406

Transparent

PP Storage Box for 1,5ml (1,8ml, 2ml) vials or 2ml shell vials, with cover (130 x 130 x 45mm), 81 cavities with alphanumeric coding of all 4 margins as well as the cavities at the bottom

Minimum order quantity: 5 pieces per colour

PP Storage Boxes




PGL190409

PP Storage Box for 1,5ml vials or 2ml shell vials, blue, with cover 67 x 67mm, 16 cavities



PGL190471

PP Storage Box for 4ml vials or 4ml shell vials, red, with cover 130 x 130 x 52mm, 49 cavities with alphanumeric coding at the cavities




PGL190648

PP Storage Box for 5ml, 10ml and 20ml Headspace-Vials, blue, with cover 130 x 130 x 102mm, 25 cavities



PGL190400

PP Vial-Rack, 200 x 105 x 17mm, for 1,5ml vials, 50 cavities, blue, stackable



PGL190472

PP Vial-Rack, 230 x 117 x 28mm, for 4ml vials, 50 cavities, blue, stackable

Minimum order quantity:
5 pieces per colour

1 pc. per PE-Bag

1 pc. per PE-Bag



HPLC and GC certified Vial KITS Certifications are getting more and more important, in order to make processes more reproducible and avoid possible sources of errors from the beginning. For La-Pha-Pack highest quality, consistency and quality control have always been very important and have been documented by HPLC and GC certified vials and closures.

- Each batch of HPLC and GC-certified KIT is tested on 15 critical parameters. In a realistic method, an HPLC/UV and GC/MS-test of vials/closure combination of blanks and contaminations will be carried out.
- The batch-specific test certificate with the HPLC and GC-Chromatograms can be handed out upon request.
- The HPLC and GC certified KITS are delivered completely shrinkwrapped for reasons of originality, purity and transport safety. This means an additional safety for the end user.

LC/MS-GC/MS certified KITS & HPLC/GC certified KITS



PGL190386

HPLC/GC Certified Vial Kit:
1.5ml Short Thread Vial, clear glass, 1st hydrolytic class, label; UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1,0mm

Pack of 100 pcs.



PGL190387

HPLC/GC Certified Vial Kit:
1.5ml Short Thread Vial, amber glass, 1st hydrolytic class, label; UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1,0mm

Pack of 100 pcs.



PGL190388

LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, clear glass, wide opening, with overwind-barrier; Ultra High Performance Seal: PP Short Thread Cap, blue, centre hole; Silicone darkblue-translucent/PTFE natural, 35° shore A, 1,0mm

Pack of 100 pcs.



PGL190389

LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, amber glass, wide opening, with overwind-barrier; Ultra High Performance Seal: PP Short Thread Cap, blue, centre hole; Silicone darkblue-translucent/PTFE natural, 35° shore A, 1,0mm

Pack of 100 pcs.

ProFill 25mm HPLC Syringe Filters



PGL190649

ProFill 25mm HPLC Syringe Filter, green, PTFE, pore size 0.2µm



PGL190650

ProFill, 25mm HPLC Syringe Filter, natural, PTFE, pore size 0.45µm



PGL190651

ProFill, 25mm HPLC Syringe Filter, blue, Regenerated Cellulose (RC), pore size 0.2µm



PGL190652

ProFill, 25mm HPLC Syringe Filter, yellow, Regenerated Cellulose (RC), pore size 0.45µm



PGL190653

ProFill, 25mm HPLC Syringe Filter, bright blue, Nylon (PA), pore size 0.2µm



PGL190654

ProFill, 25mm HPLC Syringe Filter, bright green, Nylon (PA), pore size 0.45µm

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs

Pack of 100 pcs.

Pack of 100 pcs.

17mm HPLC Syringe Filters



PGL190486

17mm HPLC Syringe Filter, PTFE, pore size 0.2µm, injection moulded ring blue + print of the membrane type



PGL190487

17mm HPLC Syringe Filter, PTFE, pore size 0.45µm, injection moulded ring yellow + print of the membrane type



PGL190488

17mm HPLC Syringe Filter, Regenerated Cellulose (RC), pore size 0.2µm, injection moulded ring granite + print of the membrane type



PGL190489

17mm HPLC Syringe Filter, Regenerated Cellulose (RC), pore size 0.45µm, injection moulded ring brown + print of the membrane type



PGL190490

17mm HPLC Syringe Filter, Nylon (PA), pore size 0.2µm, injection moulded ring purple + print of the membrane type



PGL190491

17mm HPLC Syringe Filter, Nylon (PA), pore size 0.45µm, injection moulded ring green + print of the membrane type



PGL190492

17mm HPLC Syringe Filter, Glass Fibre Prefilter/PVDF, pore size 0.2µm, injection moulded ring black + print of the membrane type



PGL190493

17mm HPLC Syringe Filter, Glass Fibre Prefilter/PVDF, pore size 0.45µm, injection moulded ring red + print of the membrane type

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.


PGL190663

30mm HPLC Syringe Filter, Glass Fibre Prefilter/PTFE, pore size 0.2µm, injection moulded ring blue + print of the membrane type

Pack of 100 pcs.


PGL190664

30mm HPLC Syringe Filter, Glass Fibre Prefilter/PTFE, pore size 0.45µm, injection moulded ring yellow + print of the membrane type

Pack of 100 pcs.


PGL190665

30mm HPLC Syringe Filter, Glass Fibre Prefilter/Regenerated Cellulose (RC), pore size 0.2µm, injection moulded ring granite + print of the membrane type

Pack of 100 pcs.


PGL190666

30mm HPLC Syringe Filter, Regenerated Cellulose (RC), pore size 0.45µm, injection moulded ring brown + print of the membrane type

Pack of 100 pcs.


PGL190667

30mm HPLC Syringe Filter, Glass Fibre Prefilter/Nylon (PA), pore size 0.2µm, injection moulded ring purple + print of the membrane type

Pack of 100 pcs.


PGL190668

30mm HPLC Syringe Filter, Glass Fibre Prefilter/Nylon (PA), pore size 0.45µm, injection moulded ring green + print of the membrane type

Pack of 100 pcs.


PGL190669

30mm HPLC Syringe Filter, Glass Fibre Prefilter/PVDF, pore size 0.2µm, injection moulded ring black + print of the membrane type

Pack of 100 pcs.


PGL190670

30mm HPLC Syringe Filter, PVDF, pore size 0.45µm, injection moulded ring red + print of the membrane type

Pack of 100 pcs.


PGL190671

30mm HPLC Syringe Filter, Glass Microfibre, pore size 1.2µm, injection moulded ring orange + print of the membrane type

Pack of 100 pcs.

Plastic Disposable Syringes with Luer Lock



PGL190009

2 ml Luer-Lok Plastic Disposable Syringe, non sterile



PGL190018

5 ml Luer-Lok Plastic Disposable Syringe, non sterile



PGL190211

10 ml Luer-Lok Plastic Disposable Syringe, non sterile



PGL190607

20 ml Luer-Lok Plastic Disposable Syringe, non sterile

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Plastic Disposable Syringes with Luer Slip



PGL190010

2 ml Luer-Slip Plastic Disposable Syringe, non sterile



PGL190019

5 ml Luer-Slip Plastic Disposable Syringe, non sterile



PGL190212

10 ml Luer-Slip Plastic Disposable Syringe, non sterile



PGL190608

20 ml Luer-Slip Plastic Disposable Syringe, non sterile

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Pack of 100 pcs.

Poznámky

pragolab



Laboratórne prístroje a vybavenie laboratórií

www.pragolab.sk