


TraceGOLD TG-5MS GC Columns

The most widely used MS phase in gas chromatography

- Low polarity phase
- Low bleed for excellent signal-to-noise ratio, sensitivity and mass spectral integrity
- Exceptional inertness ideal for analysis of active compounds
- Equivalent to USP G27 phase



TraceGOLD TG-5MS GC Columns

ID (mm)	Length (m)	Film Thickness (µm)	Cat. No.	Quantity
0.10	10	0.1	26098-0200	1 Each
0.15	10	0.15	26098-2750	1 Each
	20	0.15	26098-2760	1 Each
	40	0.15	26098-2940	1 Each
0.18	20	0.18	26098-5780	1 Each
	20 with 5m SafeGuard	0.18	26098-5785	1 Each
0.25	10	0.25	26098-1180	1 Each
	15	0.25	26098-1300	1 Each
		1.0	26098-2840	1 Each
		0.25	26098-1305	1 Each
	30	0.25	26098-1420	1 Each
		0.5	26098-2230	1 Each
		1.0	26098-2960	1 Each
		0.1	26098-0475	1 Each
	30 with 5m SafeGuard	0.25	26098-1425	1 Each
		0.25	26098-1540	1 Each
		0.5	26098-2350	1 Each
0.32	60	1.0	26098-3080	1 Each
		0.25	26098-2850	1 Each
		0.25	26098-1430	1 Each
	30	0.5	26098-2240	1 Each
		1.0	26098-2970	1 Each
		0.25	26098-1435	1 Each
0.53	15	0.25	26098-1550	1 Each
		0.5	26098-2360	1 Each
		1.0	26098-3090	1 Each
		0.25	26098-1320	1 Each
		0.5	26098-2130	1 Each
	30	1.0	26098-2860	1 Each
		1.5	26098-3340	1 Each
		0.25	26098-1440	1 Each
		0.5	26098-2250	1 Each
		1.0	26098-2980	1 Each
60	1.5	26098-3360	1 Each	
	5.0	26098-4100	1 Each	

TraceGOLD TG-5MT Metal GC Columns

ID (mm)	Length (m)	Film Thickness (µm)	Cat. No.	Quantity
0.25	15	0.1	26M98-0350	1 Each
		0.25	26M98-1300	1 Each
	30	0.25	26M98-1420	1 Each
		0.25	26M98-1540	1 Each

Applications:

- Semi-volatiles
- Phenols
- Amines
- Residual solvents and solvent impurities
- Drugs of abuse
- Pesticides
- PCB congeners
- Aroclor mixes

Similar to:

- Rxi-5ms
- DB-5
- HP-5
- HP-5ms
- Ultra-2
- SPB-5
- Equity-5
- CP-Sil 8